



## FoundationOne® CDx Specimen Guideline

Below are Specimen Guidelines to help ensure successful genomic profiling.

### Selecting the best specimen from multiple options

Has the patient been treated with targeted therapy?


<b>Yes</b>	MUST use post-targeted therapy specimen, if available.
<b>No</b>	Use recurrence or original resection (most recent material preferred). Metastasis biopsy or primary tumor acceptable (choose site with highest percent tumor or largest tumor focus).

### Acceptable Samples


- FFPE specimens, including core needle biopsies, fine-needle aspirates and effusion cytologies.
- Das Tissue should be formalin-fixed, paraffin embedded. Use standard fixation methods to preserve nucleic acid integrity. 10% neutral-buffered formalin for 6-72 hours is industry standard. DO NOT use other fixatives (Bouins, B5, AZF, Holland's).
- Do not decal. When decalcification is required, EDTA is recommended. Do not use strong acids (e.g. hydrochloric acid, sulfuric acid, picric acid).

**1 Sample Size**

When feasible, please send the block + 1 original (not recut) H&E slide.      10 unstained slides (positively charged and unbaked at 4-5 microns thick) + 1 original (not recut) H&E slide.




OR




**2 Sample Size Surface Area**

Optimal: 25 mm<sup>2</sup>      Minimal: 5 mm<sup>2</sup>

If sending slides, provide 10 unstained slides cut at 4-5 microns thick.      For small (<25mm<sup>2</sup>) or impure samples, additional unstained slides may be needed to extract sufficient DNA for testing.




OR




**3 Tumor Nuclei Percentage**

Optimal: 30%      Minimal: 20%


Prozentsatz Percent tumor nuclei = number of tumor cells divided by total number of all cells with nuclei. \*




Resection



Small Biopsy



Fine Needle Aspiration (Cell Block)



Fluid Exfoliative Cytology (Cell Block)

\* Liver specimens may require additional tumor.

